


The new label on the gako unguator jars - The new Gako quality seal



The gako quality seal offers following advantages:

- It offers a better traceability of the jar all the way back to the workers who produced each single parts and assembled them to a single jar and back to the producer of the material and color components.
- The serial number on the label enables us to react faster to complaints and to implement corrective and preventive measures.
- It's also an security feature: The combination of serial number on the jar and batch number on the packaging unit is unique, so that is impossible for criminals to falsify our jars.
- Even without the label on the packaging unit, it is immediately apparent which product it is and when the jars has been produced.
- The second date is the usability of the jar including the warranty period of 2 years. After the warranty date, the jar can be used without any loss of quality until the indicated period of use if the jars are stored at the correct conditions. Reference samples of articles, which were produced in one shift, are stored over a period of two years. This is not required for packaging material, but part of the new gako quality promise.

**Remove the label from the jar before handout the ointment to the customer.
You can use it for documentation or dispose it if not needed.**

The producer of all gako unguator products
Gako Deutschland GmbH, D-96110 Scheßlitz
is certified according to
DIN EN ISO 9001:2015

gako

Gako Deutschland GmbH
96110 Scheßlitz / Germany

Matthias Konietzko

Dipl.-Wirtsch.-Ing. Matthias Konietzko
01.02.2020